

Matlab Code

**For reading and processing
Hitachi ETG-4000 data**


A Quick Overview

*Rolf Saager, Institute of Optics,
University of Rochester*

*.m Files

- Dataprep.m
 - Executive function for data reading/processing
- Hread.m
 - Reads and parses Hitachi data
 - Also calls:
 - Probe3x3config.m
 - Probe3x5config.m
 - Probe4x4config.m
- Smoothdata.m
 - Smooths data in two ways:
 - FFT filter
 - Savitsky-Golay filter
- HbCalc.m
 - Converts Data from μ_a to Oxy-Deoxy concentration changes
 - Also calls:
 - LoadHb.m


Block Diagram of Code



Convert Raw Data to $\Delta\mu_a$

- From Beer's Law: $\frac{S}{S_o} = e^{-\Delta\mu_a^\lambda \langle L \rangle}$
- In DataPrep,
$$\Delta\mu_a^\lambda = \frac{-\log\left(\frac{S}{S_o}\right)}{\langle L \rangle}$$

- ▮ S is the timecourse
- ▮ S_o is the average value of the raw data from the end of the prescan time to the end of the initial wait period
- ▮ **Assumption:** $\langle L \rangle = 5$ for Source-Detector Separation of 3cm




Convert $\Delta\mu_a$ to HbO_2 and Hb

■ HbCalc.m

- ┌ Input
 - $\Delta\mu_a$ timecourse
 - Specific wavelength
- ┌ LoadHb.m finds the wavelength specific absorptivities, ϵ , for oxy and deoxy from a look-up table.
- ┌ Concentrations are calculated by taking the inverse of the absorptivity matrix.

$$\begin{pmatrix} \Delta\mu_a^{690nm} \\ \Delta\mu_a^{830nm} \end{pmatrix} = \begin{pmatrix} \epsilon_{HbO_2}^{690nm} & \epsilon_{Hb}^{690nm} \\ \epsilon_{HbO_2}^{830nm} & \epsilon_{Hb}^{830nm} \end{pmatrix} \begin{pmatrix} \Delta[HbO_2] \\ \Delta[Hb] \end{pmatrix}$$



Smoothing

- Two methods of smoothing, performed on each epoch
 - FFT filtering
 - 1) removes baseline frequencies of each epoch
 - 2) removes high frequency spikes (artifacts) based on a threshold
 - Savitisky-Golay
 - Smoothes data in mid frequency range
 - 2 parameters
 - Window size
 - Order of filter

Final Output

- $\text{ExpData}\{k\}.\text{HbO2}\{i\}(x,y)$
 $\text{ExpData}\{k\}.\text{Hb}\{i\}(x,y)$
where k =stim type, i =channel, x =epoch #, y =data points
- In light of the new type of stimulus
 - Data ought to be kept as an entire time course?
- Other changes:
 - Output timepoints? _____
 - Other smoothing/signal processing tools?
- Analysis tools?